



Combined Capillary Flow Porometer/ Liquid Extrusion Porosimeter



CFP-LEP-1100A

Not just products...solutions!

PRINCIPLE (Capillary Flow Porometer)

A wetting liquid is allowed to spontaneously fill the pores in the sample and a nonreacting gas is allowed to displace liquid from the pores. The gas pressure and flow rates through wet and dry samples are accurately measured. The gas pressure required to remove liquid from the pores and cause gas to flow is given by

$$D = \frac{4\gamma \cos \theta}{p}$$

where,


D is the pore diameter,

γ is the surface tension of liquid,

θ is the contact angle of liquid, and


p is the differential gas pressure.

From measured gas pressure and flow rates, the pore throat diameters, pore size distribution, and gas permeability are calculated.



DESCRIPTION (Capillary Flow Porometer)

The PMI Capillary Flow Porometer is used for R&D and quality control in industries worldwide such as filtration, nonwovens, pharmaceutical, biotechnology, healthcare, household, food, hygienic products, fuel cell, water purification, and battery. Samples often tested include filter media, membranes, paper, powders, ceramics, battery separators and health care products.



APPLICATION (Capillary Flow Porometer)

Advanced Capillary Flow Porometers yield very objective, accurate and reproducible results, considerably reduce test duration, and require minimal operator involvement. Advanced Porometers are fully automated and are designed for linear turbulence-free test gas flow. The pressure is measured close to the sample and therefore, the correction term in the differential pressure measurement is minimized. Required amount of pressure is uniformly applied on the o-ring seals on the sample and the need for hand tightening the cap on the sample chamber to apply pressure on the o-ring is eliminated. Automatic addition of wetting liquid reduces test time appreciably. This sophisticated instrument has found applications in a wide variety of industries.

TESTING CAPABILITIES (Capillary Flow Porometer)


Diameter of the most constricted part of a through pore (pore throat)

- Mean flow pore diameter (50% of flow is through pores smaller than the mean flow pore)
- Pore diameter range)
- Pore distribution:

$$f = -d[(f_w/f_d) \times 100] / dD$$



f_w = flow rate through wet sample

f_d = flow rate through dry sample



FEATURES (Capillary Flow Porometer)

- Testing of small samples as well as complete parts
- Any sample geometry (Example: sheets, rods, tubes, hollow fibers, cartridges, & powders)
- Any nonwetting liquid
- Tests in QC, research, or any number of user defined modes
- See-through sample chamber for visual observation of test available
- Real time graphic display
- Window based software for all control, measurement, data collection, data reduction, and report preparation
- Adjustable pressure on O-rings through pneumatically controlled piston-cylinder device
- Measurement of pressure close to the sample to minimize pressure drop correction
- Straight flow path avoiding turbulence
- Versatile sample chamber for a variety of samples and test modes



PMI SOFTWARE

We work closely with our customers to provide the most user friendly software for porometry. PMI Capwin software is updated to meet customer needs & requirements. The comprehensive software can be used for all PMI porometers. The software is customized to offer convenient operation with default setting for beginners & full access to all relevant measuring parameters for advanced researchers:

- Capwin manages manual instrument control, automated measuring routines ("jobs") and report print out or graph
- Capwin Data manager for interactive evaluation of measured data as well as providing sophisticated tools for creating reports & generating templates for graphs, tables & screen views
- Capwin user manager for comprehensive user management regarding user access, control & assignment of specific jobs
- On line diagnostic from anywhere in the world
- User defined definitions of paths & sub directories for data filling



SPECIFICATION (Capillary Flow Porometer)

- **Pressure Accuracy:** 0.15% of reading
- **Test Pressure:** 100, 200, and 500 psi instrument versions (700, 1400, 3500 kPa instrument versions)
- **Pressure & Flow Resolution:** 1/60,000 of full scale (1 part in 60,000)
- **Maximum Pore Size Detectable:** 500 μm
- **Minimum Pore Size Detectable:** Varies with intrusion liquid
- **Flow Rates:** Up to 200 SLPM (liters per minute)
- **Sample Geometry:** Sheets, Rods, Tubes, Hollow Fibers, Cartridges, Powders
- **Sample Sizes:**
 - Standard: 0.25" - 2.5" diameter (up to 1.5" thick)
 - Standard: 5 mm - 60 mm diameter (up to 40 mm thick)
 - Others available

| Fluid | Surface Tension dynes/cm | Diameter μm (100 psi Porometer) | Diameter μm (200 psi Porometer) | Diameter μm (500 psi Porometer) |
|----------------------|--------------------------|--|--|--|
| Water | 72 | 0.30 | 0.15 | 0.06 |
| Mineral Oil | 34.7 | 0.14 | 0.07 | 0.03 |
| Petroleum Distillate | 30 | 0.12 | 0.06 | 0.03 |
| Denatured Alcohol | 22.3 | 0.09 | 0.05 | 0.02 |
| Silwick | 20.1 | 0.08 | 0.04 | 0.02 |
| Porewick | 16 | 0.07 | 0.03 | 0.014 |
| Galwick | 15.9 | 0.07 | 0.03 | 0.014 |

DESCRIPTION (Liquid Extrusion Porosimeter)

The PMI Liquid Extrusion Porosimeter is a unique instrument with the ability to measure through-pore volume, volume distribution and liquid permeability without using mercury. The instrument is employed for characterization of porous materials used in many industries such as biotech, pharmaceutical, filtration, food, and environment. It produces no harmful effects on personnel or environment.




PRINCIPLE (Liquid Extrusion Porosimeter)

The sample is placed on a membrane in the sample chamber. The membrane is such that its largest pore is smaller than the smallest pore to be tested. The pores of the sample and the membrane are filled with a wetting liquid. The pressure of a nonreacting gas is increased on the sample to extrude the liquid from the pores. The differential pressure, p , required to displace liquid from a pore is related to its diameter, D , surface tension of the liquid, γ , and contact angle of the liquid, θ .


$$p = 4\gamma \cos \theta / D$$

The displaced liquid passes through the liquid-filled pores of the membrane and its volume is measured, while the liquid-filled pores of the membrane prevent the gas from passing through because of insufficient pressure. The gas pressure gives the pore diameter. The volume of displaced liquid gives the pore volume. Measurement of liquid flow rate without the membrane under the sample yields liquid permeability of the sample.




FEATURES (Liquid Extrusion Porosimeter)

- No toxic material like mercury is used.
- No health hazard and no disposal-related cost.
- Fully automated. Simple to use.
- Very little operator involvement.
- Highly reproducible & accurate.
- A wide variety of samples can be investigated.
- Pressure required almost an order of magnitude less than needed for mercury intrusion.
- Can be used for pressure sensitive materials.
- Only instrument capable of measuring through-pore volume.
- Capable of measuring very large pores (up to 1000 microns).
- Options includes measuring effects of application environment (stress, temperature, chemical environment) and liquid permeability.



SPECIFICATION (Liquid Extrusion Porosimeter)

- Pressure Range: 0 - 100 psi (Others Available)
- Pore Size Range: 1000 μm - 1 μm
- Resolution: 1 in 20,000
- Intrusion Volume Range: 0.01cc - 10 cc
- Sample Size: .1 to 1" diameter (.1 to 1" tall) (Others Available)



SAMPLE ANALYSIS REPORT

Liquid Extrusion Porosimetry Analysis
for
ABC Company
on 29-07-2014

by
POROUS MATERIALS, INC. ANALYTICAL SERVICES DIVISION
20 DUTCH HILL ROAD
ITHACA, NY 14850 USA
PHONE (607)-257-5544 or 1-800-TALK-PMI

File: N:\files\server1\labdata\LEP\ABC Company\Sample 2_0011-mc-c.cst
Sample ID: 2 adv (piece 2)
Operator: TLM
Lot Number: 9421
PMI Test Number: T901
Purchase Order Number: NA

Sample Weight = 4.2472 g
Bulk Density = 0.9606 g/cc
Liquid = Galwick
Liquid Extrusion Contact Angle = 0 Degrees
Liquid Extrusion Surface Tension = 16 Dynes/cm

Cumulative pore volume in cc/g of sample
% total pore volume = % of total cumulative pore volume belonging to pores of diameter > D
Average pressure = square root of $P(i) \cdot P(i-1)$
Pore size distribution function is equal to $dV/d\log P$
Surface area assumes cylindrical pores

| Pressure PSIA | Pore Diameter Microns | Cumulative Pore Vol. cc/g | %Porosity Pore Vol. PSIA | % of total Pressure $\log P$ | Average $\log P$ | dV Surf. area m ² /g | Cumulative |
|---------------|-----------------------|---------------------------|--------------------------|------------------------------|------------------|---------------------------------|------------|
| 0.007 | 1325.5128 | 0 | 0 | 0 | 0.007 | 0.051 | 0 |
| 0.008 | 1159.8078 | 0.0029 | 0.282 | 0.729 | 0.007 | 0.053 | 0 |
| 0.009 | 1030.9658 | 0.0057 | 0.545 | 1.408 | 0.008 | 0.059 | 0 |
| 0.01 | 927.7021 | 0.0084 | 0.804 | 2.079 | 0.009 | 0.074 | 0 |
| 0.014 | 656.2053 | 0.0195 | 1.872 | 4.841 | 0.012 | 0.109 | 0.0001 |
| 0.02 | 463.8572 | 0.0359 | 3.444 | 8.904 | 0.017 | 0.187 | 0.0002 |
| 0.03 | 309.2659 | 0.0688 | 6.61 | 17.09 | 0.024 | 0.335 | 0.0006 |
| 0.04 | 231.8506 | 0.1107 | 10.633 | 27.489 | 0.035 | 0.514 | 0.0012 |
| 0.05 | 185.5433 | 0.1605 | 15.422 | 39.871 | 0.045 | 0.662 | 0.0021 |
| 0.06 | 154.8203 | 0.213 | 20.46 | 52.886 | 0.055 | 0.707 | 0.0034 |
| 0.07 | 132.5505 | 0.2603 | 25.005 | 64.646 | 0.065 | 0.648 | 0.0047 |
| 0.08 | 115.9868 | 0.2979 | 29.617 | 73.982 | 0.075 | 0.557 | 0.0059 |
| 0.09 | 103.9883 | 0.3264 | 31.355 | 81.061 | 0.085 | 0.491 | 0.007 |
| 0.1 | 92.7902 | 0.3484 | 33.467 | 86.522 | 0.095 | 0.412 | 0.0079 |
| 0.11 | 84.3504 | 0.3654 | 35.103 | 90.752 | 0.105 | 0.324 | 0.0086 |
| 0.12 | 77.3277 | 0.3777 | 36.278 | 93.789 | 0.115 | 0.219 | 0.0092 |
| 0.13 | 71.3713 | 0.3853 | 37.009 | 95.679 | 0.125 | 0.152 | 0.0096 |
| 0.14 | 66.2753 | 0.3895 | 37.417 | 96.732 | 0.135 | 0.073 | 0.0099 |
| 0.15 | 61.8567 | 0.3917 | 37.626 | 97.273 | 0.145 | 0.043 | 0.01 |
| 0.2 | 46.3921 | 0.3971 | 38.144 | 98.613 | 0.173 | 0.04 | 0.0104 |
| 0.22 | 42.1755 | 0.3987 | 38.303 | 99.023 | 0.21 | 0.047 | 0.0106 |
| 0.24 | 38.6915 | 0.4005 | 38.474 | 99.455 | 0.23 | 0.062 | 0.0108 |
| 0.26 | 35.6869 | 0.4027 | 38.661 | 100 | 0.25 | 0 | 0.011 |

Summary Sheet
for
ABC Company
on 16-02-2005

File: W:\files\server1\labdata\LEP\ABC Company\Sample 2_0011-mc-c.cst
Sample ID: 2 adv (piece 2)
Operator: TLM
Lot Number: 9421
PMI Test Number: T901
Purchase Order Number: NA

Total intrusion volume = 0.4027 cc/g
Total % Porosity = 38.6807
Total Surf. area = 0.011 m²/g
Median Pore Diameter (Based on Volume) = 161.4953 Microns
Median Pore Diameter (Based on Surface Area) = 121.746 Microns
Standard Deviation (Based on Volume) = 64.8597 Microns
Standard Deviation (Based on Surface Area) = 45.6799 Microns
Average Pore Diameter (4V/S) = 146.6221 Microns
Diameter of maximum value of $dV/d\log D$ = 154.6203

Input Test Variables

Record Time Points Volume vs. Time Test Liquid Permeability Vacuum Number of Porepals

Initial Saturation Desaturation (normal test) Resaturation (hysteresis)

Sample Name: Blank
Lot Number: N/A
Run Number: 1
Operator: js

Mass of Sample: 1 g Density of Sample: 1 g/cc

File Name: C:\USERS\FPMI\DESKTOP\BLANKS\LEP-BLANK-OPEN-008.PDT

Current Group: DEFAULT
Parameter File: C:\USERS\FPMI\DESKTOP\LEP 1.5.18\DEFAULT\CSPARAM.INI

Change Group Change Parameter File

Start Test Cancel Pause at end of test

SALES & SERVICES

Our sales team is dedicated to helping our customers find which machine is right for their situation. We also offer custom machines for customers with unique needs. To find out what we can do for you, contact us. We are committed to customer support including specific service products, short response times & customer specific solutions. To quickly & flexibly meet our customer's requirement, we offer a comprehensive range of services.



CUSTOMIZE YOUR MACHINE TODAY

The most advanced, accurate, easy to use
and reproducible porometers in the world.

Celebrating



of solutions



20 Dutch Mill Rd, Ithaca, NY 14850, USA
Toll Free (US & Canada): 1-800-TALK-PMI (1-800-825-5764)
Phone: 607-257-5544 Fax: 607-257-5639

Email: info@pmiapp.com www.pmiapp.com