



The PMI Advanced
Automated Core Saturator
ACS-2K-1

Not just products...solutions!

DESCRIPTION

The PMI Automated core saturator is a completely automated system that can use various algorithms for proper saturation of multiple core samples using vacuum to remove trapped air, and pressurizing pumps to make sure the cores are completely saturated, pressure pump can control pressurizing rates precisely and maintain pressures to within a 0.2% tolerance. The Windows based program provides totally automated execution of the test, data acquisition, data storage, and data management are all performed by the user-friendly software. Operator involvement is minimal and the robust machine requires a minimal amount of care.

SPECIFICATION

- Maximum saturating pressure: 2,000 psi
- Core Length: 10"
- Core Diameter: 2.5" ID 4" OD
- Saturant: Brine, Water, oil or other liquids
- Power Requirements: 220 VAC, 50/60 Hz

SALES & SERVICES

Our sales team is dedicated to helping our customers find which machine is right for their situation. We also offer custom machines for customers with unique needs. To find out what we can do for you, contact us.

We are committed to customer support including specific service products, short response times & customer specific solutions. To quickly & flexibly meet our customer's requirement, we offer a comprehensive range of services.



Customize your machine
today!

Disclaimer : Other specifications of this product are also available. Specifications subject to change without notice. Design subject to change without notice.

The most advanced, accurate, easy to use
and reproducible Core Saturators in the world.

Celebrating



of solutions



20 Dutch Mill Rd, Ithaca, NY 14850, USA
Toll Free (US & Canada): 1-800-TALK-PMI (1-800-825-5764)
Phone: 607-257-5544 Fax: 607-257-5639

Email: info@pmiapp.com

www.pmiapp.com

